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The screenshot shows the LECO software interface with the following labeled components:

- 功能列 (Function Bar):** Located at the top, containing icons for SELECT, MULTI, SELECT ALL, MULTI FIELD, BLANK, DRIFT, CALIBRATE, RECALC, INCLUDE, EXCLUDE, DELETE, OUTPUT, IMPORT, and MORE.
- 登錄列 (Login Bar):** Located below the function bar, containing LOGIN, REPLICATES, and icons for BLANK, STANDARD, DRIFT, SAMPLE, GAS DOSE, and TEMPLATE.
- 儀表板 (Dashboard):** The central area featuring a large graph showing a peak for Oxygen (%) at 0.000274. Below the graph are several monitoring gauges: FURNACE CURRENT, COOLANT TEMPERATURE (25.0), INCOMING PRESSURE (22.0), DFC PRESSURE (mmHg), MASS FLOW (450 ccm), and FURNACE POWER.
- 功能目錄選擇列 (Function Menu):** A row of buttons for ANALYSIS, DIAGNOSTICS, SETTINGS, and INSTRUMENT.
- 測漏按鈕 (Leak Check Button):** A button labeled LEAK CHECK.
- 載氣開關按鈕 (Carrier Gas On/Off Button):** A button labeled GAS.
- 吸塵器開關 (Vacuum On/Off Button):** A button labeled VAC.
- 狀態列 (Status Bar):** A bar at the bottom showing 'Ready' and 'Analyze'.
- 分析按鈕 (Analyze Button):** A large central button labeled 'Analyze'.
- 下電極升降 (Lower Electrode Up/Down):** Buttons labeled 'DOWN' and 'UP'.
- 自動清潔按鈕 (Auto Clean Button):** A button labeled 'CLEAN'.
- 中斷按鈕 (Stop Button):** A button labeled 'STOP'.

Additional annotations include:

- 單選 (Single Select):** Points to the SELECT icon.
- 複選 (Multi Select):** Points to the MULTI icon.
- 空白校正 (Blank Calibration):** Points to the BLANK icon.
- 偏移校正 (Drift Calibration):** Points to the DRIFT icon.
- 校正 (Calibration):** Points to the CALIBRATE icon.
- 重新計算 (Recalculate):** Points to the RECALC icon.
- 語言/登錄 (Language/Login):** Points to the LOGIN icon.
- 包含 (Include):** Points to the INCLUDE icon.
- 不包含 (Exclude):** Points to the EXCLUDE icon.
- 刪除 (Delete):** Points to the DELETE icon.
- 列印/匯出 (Print/Export):** Points to the OUTPUT icon.
- 下頁 (Next Page):** Points to the MORE icon.
- 展開/摺疊 (Expand/Collapse):** Points to a circular icon on the left side of the data table.
- 卷軸 (Scrollbar):** Points to the vertical scrollbar on the right side of the data table.